



# tenda®

Bleed-No-More Paste™

## Pulmonary Support Paste

**Directions:** Administer 1-60 cc syringe 2 to 4 hours prior to event. May also administer 1-60 cc syringe 6 hours out and another 2 to 4 hours out. Minimize horse's water intake until race is complete.

**Active Ingredients:**

Honey .....41%

**Inactive Ingredients:**

Apple Cider Vinegar, Dextrose, Glycerine, Xanthan Gum, Pantothenic Acid, Yunnan Baiyao, Corn Oil.

**Precautions:** Keep out of reach of children. For equine use only. Do not mix with feed. Do not freeze. If accidentally frozen allow to thaw at room temperature for 24 hours.

Net Contents: 60 cc [2 fl. oz.]  
Item # 4001

Manufactured in  
the U.S.A by:

**Tenda Horse Products**  
Sandusky, Ohio 44870  
[www.tendahorse.com](http://www.tendahorse.com)

